

Chem 101 Laboratory

Exercise #1 Laboratory

Notebook

Measurement of Volumes and Weight: Accuracy and Precision

Name: Arfaz Hossain

Lab Section: B12

Quad: 2

Date: 23 September 2023

In-lab Notes:

Experimental Procedures:

- 1 For weighing accurately an approximate 1.5 g of NaCl, I used a weigh boat to take the NaCl and then put that (subtracting the weight of the weigh boat) on a balance. The weight is documented in the Final Documentation.
- 2 Then I transferred the NaCl to a 50.00 mL volumetric flask from the weigh boat and filled to the mark with distilled water.
- 3 I took a clean dry 50 mL beaker on a balance, documented its weight (1.1).
- 4 Then I used a 10.00 mL volumetric pipette and transferred 10.00 mL of the NaCl solution from the volumetric flask to the beaker.
- 5 The weight of the solution without the beaker can be found in the Final Documentation (1.2).

I re-did the experiment for 4 more times.

Final Measurements:

Weight of the NaCl: **1.543g**.

1.1 Weight of the Beaker is **28.501g**.

1.2 Weight of the Solution-01: **11.166g**.

2.1 Weight of the Beaker is **28.503g**.

2.2 Weight of the Solution-02: **10.156g**.

3.1 Weight of the Beaker is **28.506g**.

3.2 Weight of the Solution-03: **10.183g**.

4.1 Weight of the Beaker is **28.508g**.

4.2 Weight of the Solution-04: **10.151g**.

Lab Report

Abstract:

By measuring the volume and weight of samples, the density of a solution of NaCl was determined to be approximately **1.06 g/mL**.

Data/Results:

10.00 mL of NaCl Solution	#1	#2	#3
Weight of sample from Volumetric Pipette	11.166 g	10.156 g	10.183 g
Calculated Density	1.1166 g/mL	1.0156 g/mL	1.0183 g/mL
Average Calculated Density of NaCl Solution	1.050166667 g/mL \cong 1.06 g/mL		
Standard Deviation, σ	0.046988391001277 \cong 0.05		
%RSD for the density of NaCl Solution	4.474374638% \cong 4.47%		

Algebraic Equations:

- Density, ρ (g/mL) = $\frac{\text{Weight of NaCl solution, } m \text{ (g)}}{\text{Volume of NaCl Solution, } V \text{ (mL)}}$
- Standard Deviation, $\sigma = \sqrt{\sum \frac{(\text{Density} - \text{Average Density})^2}{(n-1)}}$
- Relative Standard Deviation, $\%RSD = \frac{\sigma}{\text{Average Density}} \times 100$

Discussion:

The density of the prepared NaCl solution was determined to be approximately **1.06 g/mL** by measuring four samples of **10.00 mL** of the NaCl solution. The %RSD calculated for these measurements was **4.47%**, indicating good precision and consistency in the measurements. The accuracy of the measurement required the use of a pipette and not a cylinder because **pipettes are designed to deliver precise volumes, ensuring that each 10.00 mL sample is consistent in volume.**

Conclusions:

In this experiment, the measured density of the NaCl solution using a volumetric pipette is **approximately 1.06 g/mL**, with a % relative standard deviation of **4.47%**. These results suggest that the measurements were consistent and reliable, and the use of a volumetric pipette allowed for accurate volume transfers.

References:

- Reimer, M. et al, Laboratory Manual, Chemistry 101, pp. 13-18. (University of Victoria: Victoria, B.C.). Fall 2022.